



Auxiliary Equipment - Drying system

1. Seamless Capsules Dryer
2. Batch Dryer
3. TM
4. Drying Tunnel

Seamless Capsules Dryer

- Seamless capsules drying device
- Improvement air-permeability by double mesh
- Discharge capsules by basket Up-Down and reverse rotation
- Capsule air-gap can be expansion drying by static pressure and air volume



Batch Dryer

- It is a good dryer for small batch size with expensive products
- Can be connected to air duct from air handling unit
- Can be installed heater for fast dry



TM

- Standard size of Capsules drying device
- Each basket can be individually controlled (Drying time, and speed)
- Basket size: Ø 650 X 838mm
- About 2,470 CFM/Basket
- Holding capacity: 25,000 caps



Drying Tunnel

- Faster Drying than open type of drying room
- No cross contamination by various capsules in same room
- Save energy for small batch size of production
- cGMP standard





Auxiliary Equipment - Drying System

1. TJS
2. Double Deck
3. Capsule Tray

TJS

- Big batch size dryer
- Can be used for continuous drying system
- Basket size: Ø 730 X 1,320 mm
- Basket lift installed
- About 3,388 CFM / Basket
- Holding capacity: 55,000 caps



Double Deck

- It is high-efficient to be installed in not enough space.
- It is easy to take in/out basket



Capsule Tray

- 758(W) X 604(D) X 35(H)mm
- Material-Polypropylene (FDA approved material, 21 CFR 177.1520)
- Loading - 1,000 ~2,000 capsules per tray
- Easy to clean / Protection from scratch and break

